

AUDIOMETRIC

SP Sonic Audiometric Booth

Model SP-05



It has allowed to reduce its weight, get a high isolation index.



The SP-05 booth has a lockable door, window great functionality.

The Audiometer booth has been designed by SP Sonic under the noise specialist using newest material in noise attenuation.

SP SONIC CO., LTD
53/8 Soi Ramindra 40, Nuanchan, Buengkum,
Bangkok 10230, Thailand Fax. 02-944-6451



Tel:02-187-1184-5

Audiometric booth

Design features

- High performance noise reduction
- Easy relocation
- Internal shelf supports audiometer and to avoid damage
- Double glazed viewing window
- Can be made with left or right hand opening
- Internal incandescent light
- Ventilation duct (option ventilation fan)
- Carpet
- Heavy duty castors
- Earth leakage circuit
- Compression door handle



This is the better of an Audiometer Booth.

Designed to allow easy access by reducing the floor, permitting a flush entrance Innovative audiometric rooms and suites for all audio logical measurements and research, combine excellent performance and versatile design.

Noise reduction of specific feature

ประสิทธิภาพการลดทอนเสียง (Noise Reduction)

Hz	125	250	500	1000	2000	4000	8000
dB	20	25	38	46	50	52	55

